




FIVALI ADJUSTABLE ANKLE STABILIZER BRACE

Instructions for Use

Please read the donning instructions and warnings carefully prior to use. To ensure that the brace works properly, follow all instructions.

CARE & WASHING INSTRUCTIONS

- Before beginning prolonged use, consult with a healthcare professional or physiotherapist to ensure you're wearing the right type of brace for your specific needs and that you're wearing it correctly.
- Ensure that the ankle brace is of the appropriate size for your foot. A brace that's too tight can impede blood circulation, while one that's too loose might not offer sufficient support.
- Periodically inspect the brace for wear and tear. Stop wearing it if it starts to lose its shape, elasticity, or support.
- Regularly clean the brace to prevent skin irritations and to maintain its effectiveness.



WARNINGS ⚠

- Neoprene supports should not be worn by individuals prone to dermatitis or allergic to neoprene products.
- If you experience pain, swelling, skin irritation or any other adverse reaction while using this product, consult your healthcare professional immediately.

SINGLE USER ONLY



1 Fasten all the straps and put your foot through the ankle brace;



2 Adjust the elasticity according to your size;



3 According to your needs to adjust the X-shaped cross straps;



4 Adjust the position of the brace plates on both sides;



5 Adjust the tightening strength of the ring buckle according to the calf circumference;



6 Final application.