




# FIVALI ADJUSTABLE ELBOW SUPPORT

## Instructions for Use



### APPLICATION INFORMATION

1. Spread the brace out flatly and place it on a flat surface with the padded side down.
2. Place your arm on the stand with your palm facing up, as shown. The brace should be located directly below the elbow.
3. Pull the Velcro straps and align it over the elbow. Fasten straps to the desired compression level. Adjust straps as needed.

### CARE & WASHING INSTRUCTIONS

- Gently hand wash the elbow brace using mild soap and warm water. Avoid harsh chemicals or abrasive cleaners that can damage the fabric or materials.
- After washing, allow the elbow brace to air dry completely. Avoid machine drying or placing it near a direct heat source, as it can cause shrinkage or deformation.
- Store the elbow brace clean and dry away from direct sunlight and moisture. It is recommended to keep it in a well-ventilated area to prevent the growth of bacteria or odors.
- Regularly inspect the elbow brace for any signs of wear and tear, such as loose threads, fraying, or damaged straps. Replace the brace if any significant damage is found to ensure proper support and protection.
- Ensure that the elbow brace fits snugly but comfortably around your elbow joint. Make necessary adjustments to the straps or fasteners to achieve a secure and proper fit. Wearing an ill-fitting brace may limit its effectiveness or cause discomfort during activity.



### WARNINGS

- Neoprene supports should not be worn by individuals prone to dermatitis or allergic to neoprene products.
- If you experience pain, swelling, skin irritation or any other adverse reaction while using this product, consult your healthcare professional immediately.

### SINGLE USER ONLY

Please read the donning instructions and warnings carefully prior to use.  
To ensure that the brace works properly, follow all instructions.