




# FIVALI SPORTS ELBOW STRAP

Instructions for Use







## **APPLICATION INFORMATION**

- 1. Position the strap through the fastener to secure the brace tothe skin;
- 2. With the fastener facing outward, place the tendon pad under the elbow, toward the outside of the forearm;
- 3. Pull the strap tight until it fits snugly.

## **CARE & WASHING INSTRUCTIONS**

- Gently hand wash the elbow brace using mild soap and warm water. Avoid harsh chemicals or abrasive cleaners that can damage the fabric or materials.
- After washing, allow the elbow brace to air dry completely. Avoid machine drying or placing it near a direct heat source, as it can cause shrinkage or deformation.
- Store the elbow brace clean and dryly away from direct sunlight and moisture. It is recommended to keep it in a well-ventilated area to prevent the growth of bacteria or odors.
- Regularly inspect the elbow brace for any signs of wear and tear, such as loose threads, fraying, or damaged straps. Replace the brace if any significant damage is found to ensure proper support and protection.
- Ensure that the elbow brace fits snugly but comfortably around your elbow joint. Make necessary adjustments to the straps or fasteners to achieve a secure and proper fit. Wearing an ill-fitting brace may limit its effectiveness or cause discomfort during activity.



## WARNINGS A

- Neoprene supports should not be worn by individuals prone to dermatitis or allergic to neoprene products.
- If you experience pain, swelling, skin irritation or any other adverse reaction while using this product, consult your healthcare professional immediately.

Please read the donning instructions and warnings carefully prior to use. To ensure that the brace works properly, follow all instructions.

### SINGLE USER ONLY