



FIVALI WRIST BRACE FOR IMMOBILIZATION

Instructions for Use

Please read the donning instructions and warnings carefully prior to use. To ensure that the brace works properly, follow all instructions.

PLEASE MAKE SURE YOUR HANDS ARE CLEAN AND DRY FOR APPLICATION



1
Loosen all straps on the wrist brace. Place your hand through the wrist brace while placing your thumb through the A as shown, adjusting to the proper position;



2
Place the splint in the position shown;



3
Thread loop straps through the D-rings;



4
Gently pull on the straps until the brace is snug against the hand. The center strap (B) serves to lock the size of the splint. NOTE: Do not tie too tightly, continue to adjust as needed for better use.



5
Final application

CARE & WASHING INSTRUCTIONS

- Make sure the brace fits snugly but comfortably around your wrist.
- Take breaks from wearing the brace throughout the day to allow your skin to breathe and prevent any potential discomfort.
- Clean your wrist brace regularly to prevent the buildup of sweat, dirt, or bacteria that can cause odors or skin irritation.
- If you have a pre-existing wrist condition or injury, consult with your doctor before purchasing and using a wrist brace to ensure it is appropriate for your specific needs.



WARNINGS ⚠

- Neoprene supports should not be worn by individuals prone to dermatitis or allergic to neoprene products.
- If you experience pain, swelling, skin irritation or any other adverse reaction while using this product, consult your healthcare professional immediately.

SINGLE USER ONLY